

Contents

1. Image handling
2. Data Handling

Code location

The python code for the app is located in src/main/python

Todo

Initializing the app

The app uses fbs to build. To load external files such as icons, csv, text files when initialising the app, an ApplicationCOntext class is used. The ctx will have attributes @cached_property that point to the external files. This is passed to the main window which itself has the application context as an attribute. Everything needs to be stored under src/main/resources/base/[folder_name].

For txt and csv files, use [folder_name]=config

For images and icons, use [folder_name]=images

InspectXrays is the main class with the following main classes as attributes

1. xray_selection_menu: *class xray_selection_menu in MenuWidgets.py*
2. widget_area_menu: *class area_menu_widget in MenuWidgets.py*
3. image_widget: *class ImageHandler in ImagingWidgets.py*

Icons:

Buttons using icons are initialised using a dictionary of icons from the app context


App context

Attributes:

1. main_window: the pyqt window to be launched by FBS
2. score_profiles: h5 files that contain the metadata used to load up the score sheet for each type of scoring method
3. image_handler_icons: dictionary of icons to be used in the image handling toolbar

InspectXrays.image_widget

Class: ImageHandler



Create method in imaging-widget
that creates a polyline Tran
control pointe and adds it to the
scene

Image Handler

→ add polyline (control-points)

(x_0, y_0)

My Scene

, → add polyline (control-points)

ImageHandler.toolbar.buttons is initialised by the ApplicationContext with passes an icon_lib to ImageHandler.toolbar

ImageHandler.toolbar.buttons is a dictionary with keys (left to right in diagram below)

1. 'Draw Polyline' :

2. 'Draw Rectangle' :

3. 'Zoom Out'

4. 'Zoom In'


5. 'Undo'

6. 'Clear Label'

7. 'Good Image Quality'

8. 'Bad Image Quality'

9. 'Annotate'



ImageHandler.annotation_options: class AnnotationModelOptions<-QWidgets


AnnotationModelOptions.score_slider_layout: class score_sliders<-QVBoxLayout

Initialising score_sliders requires a list of slider names and a list of values the slider should span.

The sliders are stored as a dictionary where the key is the title of the label next to the slider and the value is the slider itself. AnnotationModelOptions.score_slider_layout.sliders or AnnotationModelOptions.score_sliders

Selected Methods:

1. get_slider_value: calls score_sliders.get_slider_values() returns a dictionary with key slider name and value the value of the slider



AnnotationModelOptions.score_sliders['Dot Size']

AnnotationModelOptions.score_sliders['Line Width']

AnnotationModelOptions.polyline_dropdown: when a new xray is loaded by selecting a new xray from combobox_xrayid, bone annotations names corresponding to the year, xrayid and anatomical object is loaded: calls self.show_selected_annotation_bone which calls self.image_widget.image_scene.add_polyline if display_polylines_box is checked

AnnotationModelOptions.display_polylines_box: if ticked, then selecting the

AnnotationModelOptions.display_rectItem_box

AnnotationModelOptions.rectItem_dropdown: when a new xray is loaded by selecting a new xray from combobox_xrayid, bone annotations names corresponding to the year, xrayid and anatomical object is loaded: calls self.show_selected_annotation_joint which calls self.image_widget.image_scene.add_rectItem if display_rectItem_box is checked


`Self.imaging_widget.toolbar_buttons['Annotate']` is connected in `InspectXrays` to `InspectXrays.trigger_annotation_mode` which calls `IMageHandler.trigger_annotation_mode` and passes `InspectXrays.output_loc` to `ImageHandler.trigger_annotation_mode`. Hence one can pass the output)loc to the profiler.

When `self.imaging_widget.toolbar_buttons['Annotate']` is called, profilers are initialised and added to the dropdown. Clicking the dropdown triggers the profiler annotate event where one annotates in turn: this is also when the next and previous button is connected

ImageHandler.joint_annotator_panel: Class JointAnnotatorPanel

Created by clicking on self.toolbar.button['Annotate'] which calls ImageHandler.trigger_annotation_mode which initialises the output_loc of JointAnnotatorPanel, and adds it to a tab; links JointAnnotatorPanel.start_button to ImageHandler.start_annotation; links JointAnnotatorPanel.reset_button to ImageHandler.reset_annotation; also initialises an item with key being the current text of the profile_dropdown, the item is HandJointAnnotationProfiler

JointAnnotatorPanel.profile_dropdown



Click start → connects next button

Click next button \Rightarrow add a rectItem. (calls imgesave
method)

\rightarrow creates a dictionary of label coord
pair e.g

my-dict['HANDS-L-DIP3'] = [x, y, w, h] or [
 \hookrightarrow save this as.

$$\begin{bmatrix} x+w & y-h \\ x-w & y-h \\ x-w & y+h \\ x+w & y+h \end{bmatrix}$$

\rightarrow when reached end of label file, asks user to
save all annotations.

InspectXrays.xray_selection_menu

Class: xray_selection_menu

Working Directory C:\Users\amr62\PycharmProjects\aspax\ saved_data

Study ID sasuke

Xray ID R280ef790d6e4fae477648e2769bef99 (2).jfif

Date 2021 Anatomical Structure hand

+ New Study Add X-ray to Study Change save location

Scoring Method PsA-MSS

Annotation PsA-MSS

0 Measure

Joint Bone Tissue Phantom Landmark

L R N/A

Start typing bone or joint name

Save Unsure View X

Area Name	Completed	Not Completed	Unsure
0 MC1_L	0.0	1.0	0.0
1 MC2_L	0.0	1.0	0.0
2 MC3_L	0.0	1.0	0.0
3 MC4_L	0.0	1.0	0.0
4 MC5_L	0.0	1.0	0.0
5 PP1_L	0.0	1.0	0.0
6 PP2_L	0.0	1.0	0.0
7 PP3_L	0.0	1.0	0.0
8 PP4_L	0.0	1.0	0.0
9 PP5_L	0.0	1.0	0.0
10 MP2_L	0.0	1.0	0.0
11 MP3_L	0.0	1.0	0.0
12 MP4_L	0.0	1.0	0.0
13 MP5_L	0.0	1.0	0.0

InspectXrays.widget_area_menu

Class: area_menu_widget

Working Directory C:\Users\amr62\PycharmProjects\aspax\ saved_data

Study ID sasuke

Xray ID R280ef790d6e4fae477648e2769bef99 (2).jfif

Date 2021 Anatomical Structure hand

+ New Study Add X-ray to Study Change save location

Scoring Method PsA-MSS

Annotation PsA-MSS

0 Measure

Joint Bone Tissue Phantom Landmark


L R N/A

Start typing bone or joint name

Save Unsure View X

Area Name	Completed	Not Completed	Unsure
0 MC1_L	0.0	1.0	0.0
1 MC2_L	0.0	1.0	0.0
2 MC3_L	0.0	1.0	0.0
3 MC4_L	0.0	1.0	0.0
4 MC5_L	0.0	1.0	0.0
5 PP1_L	0.0	1.0	0.0
6 PP2_L	0.0	1.0	0.0
7 PP3_L	0.0	1.0	0.0
8 PP4_L	0.0	1.0	0.0
9 PP5_L	0.0	1.0	0.0
10 MP2_L	0.0	1.0	0.0
11 MP3_L	0.0	1.0	0.0
12 MP4_L	0.0	1.0	0.0
13 MP5_L	0.0	1.0	0.0

InspectXrays.xray_selection_menu class xray_selection_menu



Xray_selection_menu methods


1. **getFiles:** launches a file dialog that allows one to select a file from the system returns tuple (str path/to/filename, str filetype) also sets xray_selection_menu.temp_name to be path/to/filename
2. **getDirectory:** launches a file dialog that allows one to select a folder from the system returns str path/to/foldername also sets xray_selection_menu.temp_name to be path/to/foldername
3. **create_new_study:** calls getFiles and sets the text in xray_selection_menu.current_study_info to be the full path to the selected image file used to initialise the study:- linked to xray_selection_menu.new_study_button.clicked
4. **add_xray_to_study:** calls getFiles and sets the text in xray_selection_menu.current_file_info to be the full path to the selected image file used to initialise the study:-linked to xray_selection_menu.addXrayToStudy_button.clicked
5. **change_wd:** calls getDirectory and sets the text in xray_selection_menu.wd_info to be the full path to the selected folder used to be the output directory:-linked to xray_selection_menu.set_wdir_button.clicked

InspectXrays methods

Image Data handling methods

1. open_output_folder_selector: deprecated
2. change_wd. : reads the text in self.xray_selection_menu.wd_info and sets self.output_loc to be this, clears xray_selection_menu.combox_xrayid and xray_selection_menu.combox_studyid and calls self.display_studies:- linked to self.xray_selection_menu
3. open_study_creator : loads the image file corresponding to self.xray_selection_menu .temp_name by calling self.image_widget.load_image, sets the quality to be good by pushing the like button: creates a XrayDataCreationDialog which asks the user if the image needs to have a scoresheet.....
4. open_xray_adder. : loads the image file corresponding to self.xray_selection_menu .temp_name by calling self.image_widget.load_image, sets the quality to be good by pushing the like button: creates a XrayDataCreationDialog which asks the user if the image needs to have a scoresheet.....
5. display_xrays. : checks if there are more than 0 entries in combobox_studyid, reads the current index in the combobox_studyid and loads the metafile self.output_loc/studynname/metdata.csv, pullds out the list of filenames under metadata['file_name'] and fills combobox_xrayid with these filenames. Then sets self.xray_record to be an XrayData class initialised with the meta_loc, then calls self.load_selected_xrays
6. load_selected_xrays. : Loads up the xray whose filename is currently selected in the combobox_xrayidand displays it by calling self.image_widget.load_image. Also sets the the image quality value if this is present in the metadata file

InspectXrays.widget_area_menu area_menu_widget<-distance_menu_widget



Annotation and Score tracking

This is handled by self.update_tracking_annotation.

1. Passes [bone_name]_[side_name] to self.widget_area_menu.update_table_view along with the text on the button pressed
2. Saves the new data in the table view under self.output_loc/xray_id/annotation_tracking_[date].csv

Data Structures inside the app

Xray-handling

Xrays are handled as a XrayData class found in MenuWidgets.py. Mostly handles saving annotations to XrayData.save_loc, which is set to be the InspectXrays.output_loc. Can be initialised with an image_name,xray_id,acquisition_date or with a meta_loc

Selected methods:

1. save_bone_outline(bone_id,date,plineitem): creates a folder name self.save_loc/bones/[date] and saves the plineitem in this folder as [bone_id].txt
2. save_phantom_outline(bone_id,date,plineitem): creates a folder name self.save_loc/bones/[date] and saves the plineitem in this folder as [bone_id].txt
3. save_tissue_outline(bone_id,date,plineitem): creates a folder name self.save_loc/tissue/[date] and saves the plineitem in this folder as [bone_id].txt
4. save_patch(joint_id,date,plineitem): creates a folder name self.save_loc/joint/[date] and saves the plineitem in this folder as [bone_id].txt
5. save_landmark(landmark_id,date,plineitem): creates a folder name self.save_loc/landmark/[date] and saves the plineitem in this folder as [bone_id].txt


Selected Attributes

1. meta_table: dictionary with keys: ['acquisition_date', 'xray_id', 'organ', 'file_name', 'image_quality']


Initializing the App

Xrays are added by the initialise_xray_record method

1. Calls self.connect_sub_buttons which then calls self.display_studies, which looks up all folder names in self.output_loc and adds these to combobox_studyid, if self.output_loc is not present, then it breaks out
2. InspectXrays.output_loc is not None:



InspectXrays.save_annotation



External Files

1. **distance_menu_widget.bones_joint_list**: initialised when creating the widget. Uses the file is stored in src/main/resources/base/config/joint_list.txt
2. **xray_record.meta_table**: table containing details for all the xrays corresponding to the same study

polyline_annotate_item and bone selected : calls XrayData.save_bone_outline, target is *study_id/bone/date*

polyline_annotate_item and tissue selected: calls XrayData.save_tissue_outline, target is *study_id/tissue/date*

polyline_annotate_item and Phantom selected: calls XrayData.save_phantom_outline target is *study_id/bone/date* and file name is Phantom

Rect_annotate_item and Joint selected: calls XrayData.save_patch target is *study_id/joint/date* and file name is Phantom

polyline_annotate_item and landmark selected: vcalls XrayData.save_landmark target is *study_id/landmark/date* and file name is Phantom

InspectXrays methods

Image Data handling methods

1. open_output_folder_selector: launches
2. open_study_creator :
3. open_xray_adder:
4. display_xrays. :
5. load_selected_xrays. :

ImagingWidgehandles initializing like and dislike buttons need to reinitialise like and dislike each time an image is loaded.

Images are reloaded at the following actions

1. When creating a new study done
2. when adding a new xray to the current study done
3. When changing xray from the drop down, done

need to link this to metadata

↳ changed the initialization of xRayData and XrayStudy to have a metadata column describing image qualify

↳ need to add a handler for previous metadata.csv files that do not have an image_quality column

Changes to make:

1. Add a method to test if the metadata has an image_quality column: if not, then create one with default values of 1
2. Create a method for MainWindow that changes the image_quality status of an image in the metadata.csv files
3. Link the like, dislike button to the

Tests to run:

1. Load xrays from the aspax_small_studies to see if the handler adds the image_quality column
2. Load xrays once the image_quality handle have been chenged to see if the dislike button is set checked

Flags:

1. **MyScene.draw_poly_flag**: Bool, if true, causes mouse clicks to add points to a polyline object, exclusive with draw_rect_flag, switched by clicking on ImageHandler.toolbar.buttons[‘Draw Polyline’]
2. **MyScene.draw_rect_flag**: Bool, if true, causes mouse drag and release to draw a rectItem, exclusive with draw_rect_flag, switched by clicking on ImageHandler.toolbar.buttons[‘Draw Polyline’]

Class huo_method:

-__call__: self.extract_hand_mask: 11-12.85 s


self.get_finger_midlines: 39.59 s for Gaussian filters.

self.extract_joint_features: uses cv2.Sobel which is causing a bottleneck in memory:3.6-3.7 for each phalange, perhaps can run this once and reuse these features

self.locate_joints:


Saving a score:

InspectXrays.widget_score_menu: class score_menu_widget: initialisation requires the passing of a profile, which is loaded from the resources of the application context.
Self.score_technique : name of the score technique, appears on the tab header
Self.damage_areas : list of joints that will be score under this scoring technique: used to populate the tableView
Self.damage_types. : list types of damage used to decide the number of sliders-one slider for each damage type
Self.damage_ranges: Tuple length N with list contain min and max score values for each damage type



Loading a new score sheet

Click on xray_selection_menu.score_selector-> triggers InspectXrays.load_new_score_sheet



Looks for the profile corresponding to the currently selected item in score selector, loads the profile using load_profile

